



Feature

- Test Time
- Sample Size
- Hematocrit Range
- Operating Temp Range
- Operating Humidity Range
- Glucose Test Range
- Memory Storage
- Test Averaging

Specification

- 4 Seconds & 5 daily testing reminders
- 0.4 µL
- 30 -550%
- 5 - 45 °C
- 10 - 90%
- 1.1 - 33.3 mmol/L (20 - 600 mg/dL)
- 300 results with time and date
- General, before meal and after meal
- 7, 14, 30, 60 and 90-day average