



# CCM-Fondo



The special Ultra Low Position of the Care bed is particularly suitable for people with Alzheimer`s disease. It offers fall prevention and is useable for Homecare and hospitals.



### Full motorized bed adjustment




0-70 °

Backrest




0-30 °

Footrest




67– 640 mm  
Height adjustment



Optional  
Trendelenburg/  
Antitrendelenburg

# FCCM-FONDO

*Homecare Bed*



## As standard

- 4-section steel mesh platform base 90 x 200 cm with mattress guidance at head- and footrest
- All metal work powder coated in RAL 9023 (pearl dark grey)
- Electrically operated back rest to 70° and knee break function to 30°
- Mechanical height adjustment of footrest in 5 positions
- Electrically operated height adjustment 67-640 mm, stroke 573 mm
- 4 individual lockable castors
- cover panel light beech decor, various other decors possible on demand
- Low voltage 29 Volt motor system for extra safety
- 8 buttons handset with single fault safety and locking device
- In compliance with IEC 60601-2-52
- Protection class IPX4
- Max. Patient weight: 170 KG, mattress weight 10kg.
- Safe working load: 180 KG

## As an option

- Lifting pole with triangle hand grip
- Electrically operated Trendelenburg-/ Anti-Trendelenburg function
- Reading lamp
- Underbedlight with motion sensor for Head- and Foot end